
the all NEW RATA THE STATE OF THE STATE O

Now faster and more predictable!



Reliable, tight contacts

Drawn-wire nickel titanium plus PEEK produces consistent separating pressure longer!


Easier, more predictable contours


Marginal ridge enhancements built into the matrix bands AND rings help provide that finishing flair!

Faster, less frustrating Class IIs

The Strata-G[™] rings, bands and wedges combine for a 270° Seal[™] for ultimate flash control. The rings won't spring off the tooth, and the amazing Wide Preparation ring makes tough wide preps no problem.

270° Seal[™] for ultimate flash control!





Garrison® sectional matrix systems are the standard of care.

Tofflemire-style System

- Fails to restore proximal anatomy
- Thin contact at the marginal ridge
- Large food trap (below)
- Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease



Strata-G[™] Sectional Matrix System

- + Operator-friendly retaining system
- + Naturally contoured bands
- + Anatomically correct contacts
- + Contacts at the height of contact
- Rings produce optimal tooth separation for excellent, tight contacts





Garrison® sectional matrix systems are the standard of care.

Longer lasting by design.

Advanced materials and design combine to produce our longest lasting, most predictable separator ring ever!



SG400



SG500



SG600

Anatomy top to bottom

Five sizes of super-thin stainless steel in both dead-soft and firm ensure you have the best matrix band to restore every case.



Flexibility for a better cervical seal!









^{** +91 22 6127 7777 ** +919930006127 **} sales@nkpatelandsons.com